

# Oxygen Mask

## Medium Concentration

SINGLE - USE



MEDICAL EQUIPMENT



### Area of application

Medium concentration masks are used for administration of oxygen. The Oxygen concentration delivered is 40%-60% depending on the patient's breathing. The masks are connected directly to compressed air or oxygen supply. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. The products are for single use only

### Product range

Tube Length	Connector (to Oxygen Source)	Mask Size	Description	Reference
Medium concentration oxygen mask, no tubing				
-	-	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, No tubing	93-103Q 93-103QV (Vanilla scented mask)
		Pediatric	Pediatric Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, No tubing	93-223Q
Medium concentration oxygen mask with tubing fixed to mask				
1.8m (fixed to mask)	Standard Connector	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 1.8m (7') Oxygen safety tube fixed to mask, Standard connector	8-102Q
Medium concentration oxygen mask with tubing				
2.1m	Universal Connector	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') Oxygen safety tube, Universal connector	93-102Q
2.1m	Standard Connector	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') Oxygen safety tube, Standard connector	9-102Q 9-102QV (Vanilla scented mask)
2.1m	Universal Connector	Pediatric	Pediatric Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') Oxygen safety tube, Universal connector	93-222Q
2.1m	Standard Connector	Pediatric	Pediatric Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') Oxygen safety tube, Standard connector	9-222Q
3.0m	Standard Connector	Pediatric	Pediatric Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') Oxygen safety tube, Standard connector	9-1023QV (Vanilla scented mask)
2.1m	Standard Connector	Infant	Infant Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') Oxygen safety tube, Standard connector	9-262QV (Vanilla scented mask)

### Sterilization

Non-sterile, Shelf life: 5 years

# Oxygen Mask

## Medium Concentration

SINGLE - USE



MEDICAL EQUIPMENT

Instructions for use	
Included?	Yes

Certifications	
CE-marking	CE 1639
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485 Medical Device Single Audit Program (MDSAP)

Material specifications	
Contains latex?	No
Contains phthalates?	No
Contains Bisphenols?	No
Contains cotton?	No
Contains paper?	Packaging only
Mask	Polyvinyl chloride
Elastic strap	Polyester
Nose clip	Aluminum
Oxygen nipple connector	Polypropylene
Oxygen tube	Polyvinyl chloride
Tube connector	Polyvinyl chloride

Packaging information							
REF	Primary package (Polybag)	Secondary package/Polybag	Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Qty. pcs	Qty. pcs	Dimensions mm (LxWxD)	Volume m <sup>3</sup>		
93-103Q	1	-	50	380 x 290 x 150	0.017	2.05	5,200
93-103QV	1	-	50	380 x 290 x 150	0.017	2.05	5,200
93-223Q	1	10	50	380 x 290 x 150	0.017	3.47	5,200
8-102Q	1	10	50	380 x 290 x 310	0.034	4.0	2,400
93-102Q	1	10	50	400 x 380 x 260	0.039	4.0	2,400
9-102Q	1	10	50	400 x 380 x 260	0.039	4.0	2,400
9-102QV	1	10	50	400 x 380 x 260	0.039	4.0	2,400
93-222Q	1	10	50	400 x 380 x 260	0.039	3.28	2,400
9-222Q	1	10	50	400 x 380 x 260	0.039	3.28	2,400
9-1023QV	1	10	50	400 x 380 x 260	0.039	4.55	2,400
9-262QV	1	10	50	400 x 300 x 240	0.029	3.47	1,600

Packaging material	
Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

# Oxygen Mask

## Medium Concentration

SINGLE - USE



MEDICAL EQUIPMENT

Storage and handling	
	THIS WAY UP
	Keep away from rain
	Fragile; handle with care
	Keep away from sunlight
	Do not re-use
	Do not use if package is damaged and consult electronic instructions for use
	Consult instructions for use or use electronic instructions for use

**DISTRIBUTED BY**

**UNOQUIP GMBH**  
 Schützengraben 20  
 8200 Schaffhausen  
 Switzerland  
[info@unoquip.com](mailto:info@unoquip.com)  
[www.unoquip.com](http://www.unoquip.com)

**SPONSORED IN AUSTRALIA BY**

Unoquip Pty Ltd  
 PO Box 1117  
 Windsor VIC 3181  
 1300 55 80 79  
[info@unoquip.au](mailto:info@unoquip.au)



Ningbo Shengyurui Medical  
 Appliances Co., Ltd.  
 No. 138 Binhaisi Rd.  
 Hangzhou Bay New Zone,  
 315336, Ningbo City  
 Zhejiang Province, P.R. China



Shanghai Intl. Holding  
 Corp. GmbH  
 Eiffestrasse 80  
 20537 Hamburg  
 Germany



MedNet SWISS GmbH,  
 Bäderstrasse 18,  
 5400 Baden,  
 Switzerland